

Application story


WANNERTM
HYDRA-CELL[®] PRO
SEAL-LESS PUMP TECHNOLOGIES



Nanofiltration with Purified Water/Acetonitrile Solution

Location	Switzerland	Hydra-Cell model	G13XKSJSNETK
Type of application	Nanofiltration with Purified Water/Acetonitrile Solution	Flow rate	8-10 l/min (2.1 - 2.6 US gpm)
Liquid	Purified Water/Acetonitrile Solution	Pressure	20-30 bar (290 - 435 psi)
Application details	This customer has global manufacturing facilities to provide complex modalities such as Peptides, Lipids, Oligonucleotides, Injectables, and Small Molecules including Highly Potent. Due to the customer NDA, no further details or photographs are available.		
Advantages of Hydra-Cell pump on this application			

www.hydra-cell.co.uk

Wanner Engineering -
World Headquarters & Manufacturing
Minneapolis USA
t: (612) 332-5681
e: sales@wannereng.com

Wanner International
Hampshire UK
t: +44 (0) 1252 816847
e: sales@wannerint.com

Wanner Pumps
Shanghai CHINA
t: +86-21-6876 3700
e: sales@wannerpumps.com

Wanner Pumps
Kowloon HONG KONG
t: +852 3428 6534
e: sales@wannerpumps.com